# (XVII.) PNEUMOGRAPHY

## Anatomy of respiratory system


#### **Respiratory system**

- Airways
  - Upper airways
  - Lower airways
- Lungs

### **Respiratory muscles**

- Inspiratory muscles
  - Diaphragm
  - External intercostal muscles
- Accessory inspiratory muscles
  - Scalene and sternocleidomastoid muscles
- Expiratory muscles
  - Internal intercostal muscles; abdominal muscles

# Changes of intrapleural and intraalveolar pressure (related to atmospheric pressure) during inspiration and expiration




### Inspiration

- Active process contraction of respiratory muscles
- Decrease of intrapleural pressure
- Decrease of intraalveolar pressure
- Due to pressure gradient, air flows into lungs

### **Expiration**

- Passive process (quiet expiration) elasticity of thoracic wall and lungs
- Increase of intrapleural and intraalveolar pressure
- Air flows out of lungs

# Partial pressures of gases (mm Hg) in various parts of the respiratory system and in the circulatory system





## Chemical control of breathing

### **Chemoreceptors**

- Peripheral
- Central

- Changes of  $\mathbf{pCO}_2$  ( $\mathbf{pH}$ ) or  $\mathbf{pO}_2$
- Changes of activity of respiratory neurons





## Equipment

- two respiratory belts for registration of respiratory movements
- nose clip
- sterile mouthpiece
- Krogh respirometer
- PowerLab system

## Procedure

#### Record:

- Resting respiration (1 min)
- Respiration after a mild exertion (5 squats 10 breathing cycles)
- Respiration after an intensive exertion (30 squats 10 BC)

## **Evaluation**

Following parameters in 6 chosen breathing cycles in each recorded situation


**T**i - duration of inspiration

T<sub>e</sub> - duration of expiration

BI - duration of whole breathing cycle (breathing interval)

**Ampl** - amplitude of breathing movements

Create a table, calculate arithmetic means and standard deviations



## Statistical analysis of obtained data

Choose two sets of data which will be analyzed

Follow the procedure in textbook