

IDENTIFICATION OF AN UNKNOWN ORGANIC COMPOUND

Classification Tests

In this experiment you will attempt to identify an organic unknown from a selected group of compounds from the class of alcohols, aldehydes, ketones, carboxylic acids, phenols, or amines. You will use physical properties (mp, bp), solubility characteristics and functional group tests to identify your unknown from a list of possible compounds.

Physical Properties:

Note the physical state and color of the compound. Determine the melting point or boiling point depending on the physical state of the compound. Make sure boiling point apparatus is under the snorkel hood to prevent fumes from entering the working area.

Solubility Classification:

An organic compound is generally considered soluble if ~0.2mL of a liquid or ~0.1g of solid dissolves completely in 3mL of the appropriate solvent. Determine the solubility of your unknown in the following solvents and compare your sample to the solubility characteristics of the various classes of organic compound. Remember, these are general rules-there may be exceptions depending on other functionalities in the molecule.

Solvent	Characteristics
Soluble in Water	Lower members of homologous series of O or N containing compounds (alcohols, ketones, amines, carboxylic acids, amines having <5 C atoms)
Insoluble in water but soluble in 5% NaOH and 5% NaHCO ₃	Carboxylic acids with 5 or more carbons, phenols with electron withdrawing groups and β-diketones
Insoluble in water, insoluble in 5% NaHCO ₃ , but soluble in 5% NaOH	Phenols
Insoluble in water, insoluble in 5% NaOH, 5% NaHCO ₃ , but soluble in 5% HCl	Amines

Solubility Test Procedure:

Water: Place 0.1g of powdered solid in a small test tube and add 3 successive 0.1mL portions of water. Shake the sample after each addition of water. If the solid does not appear to dissolve, warm the test tube. If the solid dissolves, cool the sample to room temperature. If the solution remains homogeneous, the sample is considered to be soluble.


For a liquid, add 0.2mL to 3.0mL water. Shake and observe whether two phases result or a homogeneous solution. The presence of two phases indicates an insoluble sample. Do not heat.

Solubility in 5% NaOH: Proceed as above but do not heat. Note whether there is an increase in temperature or any sign of reaction. If the unknown does not appear to dissolve, it is probably a negative test. However, it may be that the sodium salt that forms has only limited solubility. If in doubt, remove some of the aqueous (without solids) and place in a small test tube. Add dilute HCl dropwise until litmus paper shows an acid solution. If any solids appear, assign the sample as “soluble in 5% NaOH”.

Solubility in 5% NaHCO₃ Proceed as for water, but do not heat. Observe whether any bubbles appear, indicating the formation of CO₂, which is characteristic of a reaction with a carboxylic acid or phenol with an EWG.

Solubility in 5% HCl: Proceed as for water, but do not heat. Note any temperature increase or any other sign of reaction. If the sample does not appear to dissolve, it is probably a negative test. If in doubt, remove some aqueous sample (without solids) and place in a small test tube. Add dilute NaOH until the sample is basic to litmus. If any solids appear assign to compound “soluble in HCl”.

Use the following chart to give a preliminary classification to your unknown



Group A

Lower members of homologous series of O or N containing compounds (alcohols, ketones, amines, carboxylic acids, amines having <5 C atoms)

Group B₁

Carboxylic acids with 5 or more carbons, phenols with electron withdrawing groups (more acidic) and β-diketones (more acidic)

Group B₂

Phenols

Group C

Basic compounds (amines)

Group D

Hydrocarbons, alkyl halides, highly alkylated (containing 5 or more carbon atoms) alcohols, aldehydes, ketones.

If you determine your compound is Group A, test the pH of the solution with pH paper. A pH < 3 is indicative of a carboxylic acid or phenol. A neutral pH is indicative of an alcohol, aldehyde, or ketone. A basic pH is indicative of an amine.

Functional Group Tests:

Test for unsaturation-Rxn with Br₂: Knowns: Add 0.2 g of the known sample (to be provided) in 2 mL of methylene chloride. Add the bromine solution dropwise until the red color remains after shaking or until 25 drops have been added. Note any color change. Place a piece of moist litmus paper (blue) in the vapor space and note any indication of HBr, which would indicate a substitution reaction instead of addition.

Test for unsaturation-Rxn with KMnO_4 : Dissolve the sample to be tested in a minimal amount of water or acetone (up to 2 mL). Add the reagent as in the previous test. The formation of a brown colored solid indicates a positive test. It is sometimes helpful to place a drop of the mixture on a piece of filter paper to more easily observe the brown color.

Test for Organic Acids: Add 5 drops or 0.1g of a known carboxylic acid (provided) to a small test tube containing ~2mL of 5% NaHCO_3 . The evolution of a gas is a positive test for an organic acid.

Test for Aldehydes and Ketones: Place approximately 1g 2,4-dinitrophenylhydrazine in a 150mL Erlenmeyer flask and add 50mL isopropanol (isopropyl alcohol, IPA). Add 1mL of the unknown and heat until the solvent boils. Remove from heat and add 1.5mL conc. HCl. Reheat to a boil for 5 minutes. Remove from heat and allow the solution to cool 5 minutes. Place the solution in a ice bath and chill for 5-10 minutes. The formation of yellow/orange to red crystals is a positive test for aldehydes or ketones. The crystals may be isolated by filtration and washed with ~5mL cold isopropyl alcohol (IPA). Suck dry and dry in an oven until dry.

Tollens Test: The Tollen's reagent is fairly unstable and must be made fresh before using. To do this, place 1 mL of 5% AgNO_3 in a clean test tube. Add 1 drop of 10% NaOH and then 2-4 drops of concentrated NH_4OH until the silver oxide dissolves. Add 1-2 drops of your unknown to your freshly prepared reagent. The formation of a silver mirror or black precipitate constitutes a positive test, which is indicative of an aldehyde. The silver will deposit so as to form a mirror only on a clean glass surface. If no reaction occurs at room temperature, warm the solution slightly in a beaker of warm water.

Fehlings Test: This test is similar to Benedict's solution in that Cu^{+2} is used to oxidize an aldehyde. Add 1mL Fehling solution A to 1mL Fehling solution B. Add 2-4 drops of sample (0.1g) and heat the mixture in hot (90-100°C) water for 5 minutes or until the formation of a precipitate or color change occurs.

Jones Oxidation: Oxidation with Dichromate This test shows a positive reaction with primary and secondary alcohols. Upon oxidation the orange dichromate is reduced to the green-blue of the Cr^{+3} ion. Tertiary alcohols do not react. To each of four test tubes, add 2 ml of 1% sodium dichromate and 5 drops of concentrated sulfuric acid and mix. Then add 10 drops of the primary alcohol to the first tube, 10 drops of the secondary alcohol to the second tube, 10 drops of the tertiary alcohol to the third tube, and 10 drops of your unknown to the fourth tube. Mix and warm the test tubes in a 40-50 deg water bath for one minute. Record any color changes.

Lucas test: The Lucas reagent is a solution of zinc chloride in concentrated hydrochloric acid. Alcohols are differentiated by the rate of reaction with the the Lucas reagent to form an alkyl chloride. Tertiary alcohols react immediately and form a cloudy solution. Secondary alcohols usually react to form a cloudy solution in five minutes, and primary alcohols do not react at room temperature. To each of four test tubes, add 5 ml of the Lucas reagent. Add 1 ml of the primary alcohol to one, 1 ml of the secondary alcohol to the next, 1 ml of the tertiary alcohol to the third, and 1 ml of your unknown to the fourth test tube. Mix with a stirring rod. If the secondary alcohol does not visibly react, heat it in a 40-50 deg water bath for 30 seconds. Record any color changes.

Iodoform Test: To each of four test tubes, add 1 ml of the primary, secondary, tertiary, and unknown alcohol respectively. Add 1 ml of water to each, and 1 ml of 10% NaOH, and 1.5 ml of the iodoform reagent (0.5 M KI/I solution). Mix and record your observations. A positive test is the disappearance of the brown color and the formation of a yellow precipitate. Both of these conditions must be met in order to consider the reaction a positive test. A positive test indicates a methyl group on the carbon bearing the carbonyl or hydroxyl group.

Test for 1° and 2° amines: For the known sample, use aniline or N-methylaniline. Place 5-8 drops of acetyl chloride in a small test tube (caution: acetyl chloride is a lachrymator (eye irritant). Use under a hood. Keep the sample container capped when transporting from the fume hood to the snorkel hood). Add 10-12 drops of the sample and observe any signs of reaction. The formation of a white solid or a white vapor trail is indicative of an amine. On completion of the test slowly pour the contents of the test tube into a beaker containing 10-20 mL of NaHCO₃ solution.

