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
# PCR (Polymerase Chain Reaction)

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## 1. Vyhledat sekvenci genu – NCBI (<http://www.ncbi.nlm.nih.gov>)

The screenshot shows the NCBI Nucleotide database search results for the gene *Mus musculus* Pou5f1. The search was performed using the keyword "Nucleotide" and returned the reference sequence NM\_013633.2. The sequence is displayed in FASTA format, and the search results are shown in a table format. The table includes fields for LOCUS, DEFINITION, ACCESSION, VERSION, KEYWORDS, SOURCE, ORGANISM, REFERENCE, AUTHORS, TITLE, and JOURNAL.

NCBI Reference Sequence: NM\_013633.2

### Mus musculus POU domain, class 5, transcription factor 1 (Pou5f1), mRNA

[Comment](#) [Features](#) [Sequence](#)

LOCUS	NM_013633	1346 bp	mRNA	linear	ROD 14-FEB-2010
DEFINITION	Mus musculus POU domain, class 5, transcription factor 1 (Pou5f1), mRNA.				
ACCESSION	NM_013633				
VERSION	NM_013633.2	GI:125490391			
KEYWORDS	.				
SOURCE	Mus musculus (house mouse)				
ORGANISM	<a href="#">Mus musculus</a> Eukaryota; Metazoa; Chordata; Craniata; Vertebrata; Euteleostomi; Mammalia; Eutheria; Euarchontoglires; Glires; Rodentia; Sciurognathi; Muroidea; Muridae; Murinae; Mus.				
REFERENCE	1 (bases 1 to 1346)				
AUTHORS	Ralston,A., Cox,B.J., Nishioka,N., Sasaki,H., Chea,E., Rugg-Gunn,P., Guo,G., Robson,P., Draper,J.S. and Rossant,J.				
TITLE	Gata3 regulates trophoblast development downstream of Tead4 and in parallel to Cdx2				
JOURNAL	Development 137 (3): 395-403 (2010)				

**Change Region Shown**

**Customize View**

**Analyze This Sequence**

- ▶ Run BLAST
- ▶ Pick Primers

**Articles about the Pou5f1 gene**

- ▶ Gata3 regulates trophoblast development downstream of Tead4 and in parallel to Cdx2. [Development.
- ▶ Stk40 links the pluripotency factor Oct4 to the Erk/MAP pathway and controls extra [Proc Natl Acad Sci U S A.
- ▶ Impairment of developmental stem cell-mediated striatal neurogenesis and pluripoter [Proc Natl Acad Sci U S A.

[RefSeq Protein Product](#)

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**PCR Template** [Reset page](#) [Save search parameters](#)

Enter accession, gi, or FASTA sequence (A refseq record is preferred) [Clear](#)

NM\_013633.2

Or, upload FASTA file

**Range**

Forward primer  From  To

Reverse primer

**Primer Parameters**

Use my own forward primer (5'→3' on plus strand)  [Clear](#)

Use my own reverse primer (5'→3' on minus strand)  [Clear](#)

PCR product size

Min  Max

# of primers to return

Primer melting temperatures (T<sub>m</sub>)

Min  Opt  Max  Max T<sub>m</sub> difference

Please note the recent change in default T<sub>m</sub> calculation [?](#)

**Exon/intron selection** A refseq mRNA sequence as PCR template input is required for options in the section [?](#)

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



Total Na[+] Equivalent: 176.5

Excessive difference between product and primer melting temperatures.  
3'-end Upper Primer dimer.

*A pak už je to hračka*