Practical NMR Spectroscopy of Biomolecules

Protocols from the practical sessions have to be submitted at least 2 days before the examination!

Questions for Examination

A - Theory

- 1. Requirements for NMR samples volume, concentration, purity, solvents, buffers, salt, ions, locking and reference substances; types of NMR sample tubes
- 2. Why is magnet homogeneity important in NMR spectroscopy and how is it achieved? Methods of magnet shimming.
- 3. What is the purpose of the field/frequency lock and how does it work? How do you set and optimize the lock parameters?
- 4. Pulse calibration in direct and indirect channels. Use of the prosol table. Recalculation of the power for different pulse lengths and shapes.
- 5. Sample temperature control and calibration. What are suitable calibration samples? Explain the calibration procedure.
- 6. Describe the set-up of proton 1D measurement. What parameters should be set? Explain the interdependence among spectral width, number of points and acquisition time (parameters SW, TD, AQ).
- 7. Describe the set-up of one-dimensional ¹³C spectrum. What is the purpose of decoupling, how it works and what parameters have to be set?
- 8. Name the basic water suppression techniques and explain their principles.
- 9. Acquisition schemes in 2D spectroscopy States, TPPI, Echo-Antiecho. Explain the purpose and differences.
- 10. Homonuclear through-bond 2D correlation experiments (COSY, TOCSY). Explain the principles and applications.
- 11. Homonuclear through-space 2D correlation experiments (NOESY, ROESY). Explain the principle and applicability.
- 12. Heteronuclear 2D correlation experiments (HMQC, HSQC). Explain the principles, differences and applications.
- 13. Heteronuclear 2D correlation experiments for isotopically labeled samples. Explain the purpose and principles of sensitivity enhancement by preservation of equivalent pathways and of constant time t1 evolution.
- 14. Protein NMR spectra. Identify the signal ranges in proton, ¹³C and ¹⁵N dimensions (amide, alpha, beta, sidechain, and carbonyl).
- 15. Proton NMR spectra of nucleic acids. Identify the regions of methyl, H1', other sugar protons, pyrimidine H5, other base protons, amino, and imino signals.

- B practical measurement
 - 1. Measure proton 1D spectrum of an organic compound (0.1% Ethylbenzene in CDCl₃ standard sample). Concentrate on achieving high resolution and sensitivity.
 - 2. Measure ¹³C spectrum of an organic compound (20% Ethylbenzene in CDCl₃ standard sample) with proton decoupling.
 - 3. Measure ³¹P spectrum with proton decoupling of a DNA sample in phosphate buffer. Verify the ³¹P 90⁰ pulse length on the buffer signal.
 - 4. Measure proton 1D spectrum of the standard sucrose sample using presaturation. Evaluate the signal-to-noise ratio and resolution of the resulting spectrum on the signal of anomeric proton (5.4 ppm), use macro suppcal.
 - 5. Measure proton 1D spectrum of a DNA sample in 90% $H_2O/10\%$ D₂O. Use a method suitable for detecting imino signals.
 - 6. Set-up the 2D COSY experiment on a sample of DNA in D_2O . Verify the setup by measuring the first increment. Perform the processing, including the phase correction, on a supplied data set measured previously.
 - 7. Set-up the 2D TOCSY experiment on a sample of DNA in D_2O . Verify the setup by measuring the first increment. Perform the processing, including the phase correction, on a supplied data set measured previously.
 - 8. Set-up the 2D NOESY experiment on a sample of DNA in D₂O. Verify the setup by measuring the first increment. Perform the processing, including the phase correction, on a supplied data set measured previously.
 - 9. Measure and process ¹H-¹⁵N correlation spectrum of the amide region of ¹⁵N and ¹³C labeled Ubiquitin sample.
 - 10. Measure and process ¹H-¹³C correlation spectrum of the aliphatic region of ¹⁵N and ¹³C labeled Ubiquitin sample.