



# Quantitation Report

## Experiment Information

Run Name	Run 2021-12-15-sk1_ACO3
Run Start	15.12.2021 10:53:27
Run Finish	15.12.2021 12:30:12
Operator	MP
Notes	
Run On Software Version	Rotor-Gene Q Software 2.0.2.4
Run Signature	The Run Signature is valid.
Gain Green	10,
Machine Serial No.	0911113

## Quantitation Information


Threshold	0,030
Left Threshold	5,000
Standard Curve Imported	No
Standard Curve (1)	conc= 10 <sup>^</sup> (-0,513*CT + 13,442)
Standard Curve (2)	CT = -1,950*log(conc) + 26,215
Reaction efficiency (*)	(* = 10 <sup>^</sup> (-1/m) - 1) 2,25638
M	-1,95033
B	26,21531
R Value	0,96098
R <sup>^</sup> 2 Value	0,92348
Start normalising from cycle	1
Noise Slope Correction	Yes
No Template Control Threshold	% 0
Reaction Efficiency Threshold	Disabled
Normalisation Method	Dynamic Tube Normalisation
Digital Filter	Light
Sample Page	Page 2
Imported Analysis Settings	

## Profile


Cycle	Cycle Point
Hold 1	Hold @ 95°C, 10min 0s
Cycling (35 repeats)	Step 1: Hold @ 95°C, 15s
	Step 2: Hold @ 60°C, 30s, acquiring to Cycling A([Green][1][1])
	Step 3: Hold @ 72°C, 20s

Hold 2	Hold @ 72°C, 1min 0s
Melt	Ramp from 55°C to 90°C
	Hold for 90s on the 1st step
	Hold for 5s on next steps,Melt A([Green][1][1])


### Raw Data For Cycling A.Green



### Quantitation data for Cycling A.Green



### Standard Curve



No.	Color	Name	Type	Ct	Ct Comment	Given Conc (IU/ml)	Calc Conc (IU/ml)
21		T6 2 10X	Standard	21,17		1 000	385
22		T6 2 10X	Standard	21,23		1 000	359
23		T6 2 1x	Standard	18,70		10 000	7 170
24		T6 2 1x	Standard	18,36		10 000	10 710
25		T6 2 100x	Standard	20,93		100	515
26		T6 2 100x	Standard	21,06		100	442
27		T6 2 1000x	Standard	23,89		10	16
28		T6 2 1000x	Standard	24,01		10	13
29		T6 2 10000x	Standard	27,18		1	
30		T6 2 10000x	Standard	26,63		1	1
31		2 NK	Unknown	27,17			
32		2 NK	Unknown	26,32			1
33		D6 2 ACO2	Unknown	32,81			
34		D6 2 ACO2	Unknown	29,12			
35		D6 2 ACO2	Unknown	24,55			7
36		2 NK	Unknown	27,16			
37		B6 2 ACO2	Unknown	14,58			925 002
38		B6 2 ACO2	Unknown	14,58			929 261
39		B6 2 ACO2	Unknown	14,43			1 101 722
40		2 NK	Unknown	26,12			1

**Warning: The following samples were not analysed :**

1T6 1 1X- 2T6 1 1X- 3T6 1 10x- 4T6 1 10x- 5T6 1 100x- 6T6 1 100x- 7T6 1 1000x- 8T6 1 1000x- 9T6 1 10000x- 10T6 1 10000x- 111 NK- 121 NK- 13D6 1 ACO3- 14D6 1 ACO3- 15D6 1 ACO3- 16B6 1 ACO3- 17B6 1 ACO3- 18B6 1 ACO3- 191 NK- 201 NK- 41T6 3 1X- 42T6 3 1X- 43T6 3 10x- 44T6 3 10x- 45T6 3 100x- 46T6 3 100x- 47T6 3 1000x- 48T6 3 1000x- 49T6 3 10000x- 50T6 3 10000x- 513 NK- 523 NK- 53D6 3 ACO4- 54D6 3 ACO4- 55D6 3 ACO4- 563 NK- 57B6 3 ACO4- 58B6 3 ACO4- 59B6 3 ACO4- 603 NK

**Legend:**

NEG (NTC) - Sample cancelled due to NTC Threshold.

NEG (R. Eff) - Sample cancelled as efficiency less than reaction efficiency threshold.

