

**Topic P10: Basics of clinical mycology**

To study: Fungi

From spring term: Microscopy, culture, antibiotic susceptibility, precipitation

**Table for major results of Task 1 to Task 3 (to be filled step by step):**

Strain		K	L	M	N
Gram stain – Task 1					
Culture (blood agar) Task 2a	Size				
	Colour				
	Shape				
	Profile				
	Haemolysis				
	Surface				
	Odour				
Task 2b: growth on Sabouraud + chloramphenicol (growth +/-)					
Task 2c: growth on a chromogenic medium					
<b>FINAL CONCLUSION (according to Task 2c/Task 3)</b>					

**Task 1: Microscopy of strains of bacteria and yeasts**

Gram stain given cultures of microorganisms. Use immersion microscopy (immersion objective magnifying 100x). Write down your results to the table. Remark size differences between yeasts and bacteria

**Task 2: Culture of bacteria and yeasts**

**a) Culture on blood agar**

Describe colonies of given strains on blood agar and fill in the main table. Do not forget to describe odour. Remark, that colonies of yeasts (according to Task 1) are similar to colonies of some bacteria (especially to the colonies of G+ cocci, likely to be staphylococci according to the morphology of colonies).

**b) Culture on Sabouraud agar with chloramphenicol**

Evaluate growth of given strains on selective agar for yeasts and molds (Sabouraud agar with chloraphenicol). Sabouraud agar itself is not selective, but it is made selective with help of a broad-spectre antibiotic (chloramphenicol).

**c) Culture on a chromogenic medium**

Chromogenic media for yeasts enable differentiation of the most important species of *Candida* genus. With the help of control strains, try to assess species of *Candida* using the chromogenic medium. If the colonies are white (no colour substance present), it means, that the strain cannot be differenciated using this chromogenic medium.

**Task 3: Auxacolor**

Read the results of Auxacolor for the strain that could not be differenciated using chromogenic medium. In this case, we do not count the code, but we compare the results with the table. (You might get more than one taxon. If so, remember, that *C. albicans* would be green in 2c) and *Rhodotorula* would be pigmented (red) in 2b.)

C Neg*	Glu**	Mal**	Sac**	Gal**	Lac**	Raf**	Ino**
Cel**	Tre**	Ado**	Mel**	Xyl**	Ara**	Act***	Pox#
Result of identification:							

\*normally blue

\*\*\*yellow positive, colourless negative

\*\*yellow positive, blue negative

#brown positive, colourless negative

**Task 4: Assessment of antimicrobial drugs susceptibility**

For treatment of fungal infections, it is not possible to use antibiotics. We have to use special drugs – antimycotics. These, on the other hand, are not effective against bacteria.

**a) Assessment of susceptibility to antimycotics**

Perform in vitro susceptibility testing of given strains to antimycotics. Into the table, write the full name of the antimycotics according to a card and for all tested strains assess susceptibility or resistance. For Amphotericin B, the reference zone is 10 mm. For other antimycotics, it is 20 mm, but is not necessary, that the inside of the zone is absolutely clear.

Strain →				
Antimycotic (full name)	Zone Ø (mm)	Interpretation	Zone Ø (mm)	Interpretation

Strains \_\_\_ and \_\_\_ (i. e. strains of \_\_\_\_\_) are resistant to all given drugs.

**b) Assessment of susceptibility to antibiotics**




Perform in vitro susceptibility testing of given strains to antibiotics. Into the table, write the abbreviation of the antibiotics according to a card and for all tested strains measure the susceptibility zones. On your card, you have limit zones – according to them, interpret the zones as susceptible (S) resistant (R) and dubious (D).

Strain →				
Antibiotic (full name)	Zone Ø (mm)	Interpretation	Zone Ø (mm)	Interpretation

Strains \_\_\_ and \_\_\_ (i. e. strains of \_\_\_\_\_) are resistant to all given drugs.

**Task 5: Microscopy of molds**

Molds are usually microscopied differently than yeasts. Gram staining is used rarely. **Mostly we use wet mount, objectives magnifying 10× to 40×.** Draw and describe three given species of molds. Do not forget to describe hyphae, macro- and microconidia and other observed objects. Use lines or numbers to connect description to the objects.



Genus _____ Species _____  Description: _____ _____ _____	Genus _____ Species _____  Description: _____ _____ _____	Genus _____ Species _____  Description: _____ _____ _____
-------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------

**Task 6: Culture of molds**

Molds usually require longer time to grow. That is why we mostly do not use Petri dishes for culture, but only test tubes (to avoid drying and contamination. Draw culture results of given molds.

**Task 7: Indirect diagnostics of aspergilosis**

Evaluate results of precipitation in gel in order to diagnose antibodies in aspergilosis. Draw the result.



**Task 8: Sampling for mycoses**

Look at the video „Sampling of cutaneous mycoses. Write down basic rules for sampling in superficial mycoses

**Check-up questions:**

1. Name some antimycotics that could be used for surface mycoses
  
  
  
  
  
  
  
  
  
  
2. Name some antimycotis suitable for generalized candidoses and other systemic mycoses
  
  
  
  
  
  
  
  
  
  
3. What is the growth time (in days, approx.) of
  - a) yeasts
  - b) zygomycetes
  - c) aspergilli
  - d) dermatophytes
  
  
  
  
  
  
  
  
  
  
4. What are suitable culture media and temperatures for pathogenic fungi?
  
  
  
  
  
  
  
  
  
  
5. Name at least five factors helping systemic mycoses to get formed: