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



## Colony Morphology: Describing Bacterial Colonies

Frequently during the semester you will need to describe bacterial (or fungal) growth observed on slants or Petri plates. It will be useful to learn the terminology used for describing common colony types. The following outline will be helpful for verbally communicating the appearance of observed colonial growth.


**1. Form** – The form refers to the shape of the colony. These forms represent the most common colony shapes you are likely to encounter.


CIRCULAR

IRREGULAR

FILAMENTOUS

RHIZOID


**1a. Size** – The size of the colony can be a useful characteristic for identification. The diameter of a representative colony may be measured. Tiny colonies are referred to as **punctiform**.

**1b. Surface** – Bacterial colonies are frequently shiny and smooth in appearance. Other surface descriptions might be: **veined, rough, dull, wrinkled (or shriveled), glistening**.

**1c. Texture** – Several terms that may be appropriate for describing the texture or consistency of bacterial growth are: **dry, moist, mucoid, brittle, viscous, butyrous (buttery)**.

**1d. Color** – It is important to describe the color or pigment of the colony. Also include descriptive terms for any other relevant optical characteristics such as: **opaque, cloudy, translucent, iridescent**.

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Bacillus mycoides filamentous colony

**2. Elevation** – This describes the “side view” of a colony. These are the most common.

**FLAT**

**RAISED**

**UMBONATE**  
(having a knobby protuberance)

**CRATERIFORM**

**CONVEX**

**PULVINATE**  
(cushion-shaped)

**3. Margin** – The margin or edge of a colony (or any growth) may be an important characteristic in identifying an organisms. Several examples are shown below.

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
*Streptococcus pyogenes*(1)

kvasinky

*Candida Albicans* is a type of yeast that can grow on the surface of skin(7)

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










# Ziehl–Neelsen stain (acid fast staining)

The Ziehl–Neelsen stain, also known as the acid-fast stain, was first described by two German doctors: the bacteriologist Franz Ziehl (1859–1926) and the pathologist Friedrich Neelsen (1854–1898).

It is a special bacteriological stain used to identify acid-fast organisms, mainly Mycobacteria. Mycobacterium tuberculosis is the most important of this group because it is responsible for tuberculosis .

Acid fast organisms like Mycobacterium contain large amounts of lipid substances within their cell walls called **mycolic acids**. These acids resist staining by ordinary methods such as a Gram stain. It can also be used to stain a few other bacteria, such as Nocardia. The reagents used are Ziehl–Neelsen **carbol fuchsin, acid alcohol, and methylene blue**. Acid-fast bacilli will be bright red after staining.

**Initially, Carbol Fuchsin stains every cell. When they are destained with acid-alcohol, only non-acid-fast bacteria get destained since they do not have a thick, waxy lipid layer like acid-fast bacteria. When counter stain is applied, non-acid-fast bacteria pick it up and become blue when viewed under the microscope. Acid-fast bacteria retain Carbol Fuchsin so they appear red.**

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
Ashby's Mannitol Agar is used for the cultivation of the Azotobacter species that can use mannitol and atmospheric nitrogen as a source of carbon and nitrogen, respectively. Besides having the ability to fix atmospheric nitrogen, Azotobacter can also synthesize biologically active substances which improve seed germination and plant growth.

**Selective media** allow the growth of certain type of organisms, while inhibiting the growth of other organisms. Selective inhibition of some types of microorganisms can be studied by adding certain **dyes, antibiotics, salts or specific inhibitors** that will affect the metabolism or enzymatic systems of the organisms. For example, media containing **potassium tellurite, sodium azide or thallium acetate** at different concentrations of 0.1 - 0.5 g/l will inhibit the growth of all Gram-negative bacteria. Media supplemented with the antibiotic **penicillin** concentration 5-50 units/ml or **crystal violet** 2 mg/l inhibit the growth of Gram-positive bacteria. Tellurite agar, is used to select for Gram-positive organisms, and nutrient agar supplemented with the antibiotic penicillin can be used to select for the growth of Gram negative organisms.

**Differential media** are widely used for differentiating closely related organisms or groups of organisms. Because of the presence of certain dyes or chemicals in the media, the organisms will produce certain **characteristic changes or growth patterns that are used for identification or differentiation of microorganism.**

**Enriched media** are media that have been supplemented with highly nutritious materials such as blood, serum or yeast extract for the purpose of cultivating fastidious organisms.

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

Media	Classification	Selective and Differential agents	Type of organisms isolated
<b>Mannitol salt agar (MSA)</b>	Selective and differential	7.5% NaCl and mannitol for isolation and identification of most <i>S.aureus</i> strains	Staphylococci and Micrococci
<b>MacConkey's Agar</b>	Selective and differential	Lactose, bile salts, neutral red, and crystal violet	Gram-negative enteric bacilli
<b>Eosin methylene blue agar (EMB)</b>	Selective and differential	Lactose, eosin Y, and methylene blue	Gram-negative enteric bacilli
<b>Phenylethyl Alcohol Agar (PEA)</b>	Selective	Phenylethyl alcohol (inhibits gram negatives)	Gram-positive bacteria
<b>Hektoen Enteric Agar (HE)</b>	Selective and differential	Lactose, sucrose, bile salts, ferric ammonium sulfate, sodium thiosulfate, bromthymol blue, acid fuchsin	Salmonella and Shigella species (enteric pathogens)
<b>Blood Agar</b>	Enriched and differential	5% defibrinated sheep blood	Almost all bacteria; differential for hemolytic organisms
<b>Chocolate Agar</b>	Enriched	1% hemoglobin and supplements	Most fastidious pathogens such as Neisseria and Haemophilus

## Mannitol salt agar (MSA)

- Mannitol salt agar is both a selective and differential media used for the isolation of pathogenic *Staphylococci* from mixed cultures.
- On MSA, only pathogenic *Staphylococcus aureus* produces small colonies surrounded by yellow zones. The reason for this color change is that *S. aureus* have the ability to ferment the mannitol, producing an acid, which, in turn, changes the indicator color from red to yellow. The growth of other types of bacteria is usually inhibited. This growth differentiates *S.aureus* from *S.epidermidis*, which forms colonies with red zones or both zones.

Figure1 : Mannitol Salt Agar

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
# MacConkey's agar (MAC)

- MacConkey's Agar is both a selective and differential media; it is **selective for Gram negative bacteria and can differentiate those bacteria that have the ability to ferment lactose**
- Neutral red pH indicator - Stains microbes fermenting lactose (hot pink in acid pH, rose in neutral pH, brown in alkaline pH)
- By **utilizing the available lactose** in the medium, Lac<sup>+</sup> (lactose fermenting) bacteria such as *Escherichia coli*, *Enterobacter* and *Klebsiella* will **produce acid** in the medium, which lowers the pH of the agar below 6.8 and results in the appearance of **red or pink colonies**. Non-lactose fermenting bacteria such as, *Proteus species*, *Salmonella*, *Pseudomonas aeruginosa* and *Shigella* cannot utilize lactose in the medium, and will use peptone instead. This results in the formation of ammonia, which raises the pH of the agar, and leads to the formation of white or colorless colonies in the plate.

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
## Eosin methylene blue agar (EMB), Endo agar


- Eosin methylene blue agar (EMB) is both a selective and differential medium used for the **detection and isolation of Gram-negative intestinal pathogens**.
- **Acid production from lactose fermentation** causes precipitation of the dyes on the surface of the colony resulting in different colors.
- **Large amounts of acid - green metallic sheen**
- **Small amounts of acid - pink**
- **No fermentation - colorless**
- Enterobacter aerogenes produces large colonies which are pink-to-buff around dark centers. Escherichia coli produce small, dark colonies with a green metallic sheen. Pseudomonas, Proteus, Salmonella and Shigella sp produces colorless colonies because it does not ferment lactose.

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
# Phenylethyl alcohol agar (PEA)

- Phenylethyl Alcohol (PEA) Agar with or without 5% sheep blood is a selective medium for the **isolation of gram-positive organisms**, particularly gram-positive cocci, from specimens of mixed gram-positive and gram-negative flora.
- **Phenylethyl alcohol** – Inhibits the growth of Gram negatives since it selectively and reversibly inhibits DNA synthesis, thus selecting for Gram positives.

PEA alters the membrane permeability, of Gram-negative bacteria allowing influx of otherwise blocked molecules. This results in leakage of large amounts of cellular potassium that ultimately results in disruption or inhibition of DNA synthesis of Gram-negative bacteria.

# Hektoen Enteric Agar (HE)

- Hektoen Enteric (HE) Agar is a moderately selective medium used in qualitative procedures for the isolation and cultivation of **gram-negative enteric microorganisms**, especially *Shigella* and *Salmonella* from a variety of clinical and nonclinical specimens
- **Coliforms** capable of overcoming the moderately inhibitory qualities of the media will develop into **orange or salmon-pink colonies** in the presence of the bromthymol blue indicator. *Shigella* species develop into **green-colored colonies** with darker blue-green centers. *Salmonella* species appear as **blue-green colonies** with or without black centers. **Producers of H<sub>2</sub>S** will form **black-centered colonies** in the presence of the ferric ammonium citrate indicator.

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
## Blood agar

- Blood agar is both **differential and enriched** medium. The blood that is incorporated into this medium is an enrichment ingredient for the cultivation of fastidious organisms such as the *Streptococcus* species.
- A number of streptococcal species produce substances that destroy red blood cells; that is, they cause lysis of the red cell wall with subsequent release of hemoglobin. Such substances are referred to as **hemolysins**. The activity of streptococcal hemolysins also known as streptolysins can be readily observed when the organisms are growing on a blood agar plate.

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## Chocolate agar

- Fastidious organisms such as *Haemophilus* and *Neisseria* require specially enriched culture media and microaerophilic incubation conditions. “Chocolate” agar is commonly used for primary isolation of *Haemophilus* from clinical specimens. This medium contains **hemoglobin** derived from bovine red blood cells as well as other enrichment growth factors. Chocolate agar may be made selective for *Haemophilus* species by the addition of bacitracin.
- Two special growth factors, called X and V, are required by some *Haemophilus* species. **The X factor is hemin**, a heat-stable derivative of hemoglobin. **The red blood cells in chocolate agar have been heated until they are lysed, producing the characteristic brown color of this medium.**

Figure 7 :Uninoculated Chocolate Agar Plate

Figure 22:*Haemophilus influenzae*

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
Carbohydrate Fermentation Broth  
trypticase, Sodium chloride, and Phenol red

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

Medium: Sodium phosphate buffer,  
1 M, pH 7.0

O-Nitrophenyl- $\beta$ -D-  
galactopyranoside (ONPG), 0.75 M  
Physiologic saline

Toulene

O-Nitrophenyl- $\beta$ -D-galactopyranoside (ONPG) is structurally similar to lactose (i.e. ONPG is an analog of lactose), except that orthonitrophenyl has been substituted for glucose.

**Lactose fermenter (ONPG Positive):** *E.coli*, *Klebsiella spp*, *Enterobacter spp* produce  $\beta$ -galactosidase

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
### MRVP broth (pH 6.9)

Ingredients per liter of deionized water:  
buffered peptone= 7.0 gm  
glucose= 5.0 gm  
dipotassium phosphate= 5.0 gm

Some bacteria have the ability to utilize glucose and convert it to a stable acid like lactic acid, acetic acid or formic acid as the end product.

These bacteria initially metabolise **glucose to pyruvic acid**, which is further metabolized through the '**mixed acid pathway** to produce the **stable acid**. The type of acid produced differs from species to species and depends on the specific enzymatic pathways present in the bacteria. The acid so produced decreases the pH to 4.5 or below, which is indicated by a change in the colour of methyl red from **yellow to red**.

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

glucose= 5.0 gm

dipotassium phosphate= 5.0 gm

The Voges-Proskauer (VP) test is used to determine if an organism produces **acetylmethyl carbinol** from glucose fermentation. If present, **acetylmethyl carbinol** is converted to **diacetyl** in the presence of  $\alpha$ -**naphthol**, strong alkali (**40% KOH**), and atmospheric oxygen. The  $\alpha$ -**naphthol** was not part of the original procedure but was found to act as a color intensifier by Barritt and must be added first.

The **diacetyl** and **guanidine**-containing compounds found in the **peptones** of the broth then condense to form a **pinkish red polymer**.

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

Citrate = oxaloacetate + acetate

oxalacetate = pyruvate + CO<sub>2</sub>


**A. Under alkaline conditions**, pyruvate is metabolized to acetate and formate.

pyruvate = acetate + formate

**B. At pH 7.0 and below**, lactate and acetoin are also produced.

pyruvate = acetate + lactate + CO<sub>2</sub>

pyruvate = acetoin + CO<sub>2</sub>

Citrate utilization test

Further metabolic breakdown is **dependent upon the pH of the medium**.

**A. Under alkaline conditions**, pyruvate is metabolized to acetate and formate.

pyruvate = acetate + formate

**B. At pH 7.0 and below**, lactate and acetoin are also produced.

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
Organisms which produce the enzyme thiosulfate reductase can reduce sulfur to hydrogen sulfide gas. This happens when the strain either degrades the amino acid cysteine during protein degradation, or when anaerobic respiration shuttles the electrons to sulfur instead of to oxygen. In the SIM tubes, the medium contains casein and animal proteins as amino acid sources, sodium thiosulfate as a source of sulfur and ferrous ammonium sulfate as the H<sub>2</sub>S indicator. Cysteine is a sulfur containing amino acid present in the SIM medium.

The enzymes cysteine desulfurase and thiosulfate reductase catalyze hydrolysis reactions that produce H<sub>2</sub>S. This gas combines with the ferrous ammonium sulfate forming an insoluble, black ferrous sulfide precipitate. The black color thus acts as an indicator for the presence of hydrogen sulfide.

In the laboratory, a fresh culture of the organism is inoculated with a single stab using straight needle through the center of the medium. Following incubation, the tube is observed for H<sub>2</sub>S production (blackening of the medium).

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

**Tryptophan** is an amino acid that can undergo deamination and hydrolysis by bacteria that express tryptophanase enzyme. **Indole** is generated by reductive **deamination** from **tryptophan** via the intermediate molecule **indolepyruvic acid**. **Tryptophanase** catalyzes the deamination reaction, during which the **amine (-NH<sub>2</sub>)** group of the **tryptophan** molecule is removed. Final products of the reaction are **indole**, **pyruvic acid**, **ammonium (NH<sub>4</sub><sup>+</sup>)** and **energy**. **Pyridoxal phosphate** is required as a coenzyme.

When **indole** is combined with **Kovac's Reagent** (which contains hydrochloric acid and p-dimethylaminobenzaldehyde in amyl alcohol) the solution turns from yellow to **cherry red**. Because amyl alcohol is not water soluble, the red coloration will form in an **oily layer at the top of the broth**.

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
Phenylalanine deaminase test also known as phenylpyruvic acid (PPA) test is used to test the ability of an organism to produce enzyme deaminase. This enzyme removes the amine group from the amino acid phenylalanine and produces phenylpyruvic acid (PPA) and ammonia i.e. oxidative deamination of phenylalanine. Phenylpyruvic acid reacts with ferric iron (10% ferric chloride is added in the medium) producing a visible green color.

## Results

- 1. Positive test:** Production of green colour (Phenylpyruvic acid thus formed reacts with ferric chloride producing a green colored compound thus turning the medium dark green). *Proteus sp.*, *Morganella sp.*, *Providencia sp* give positive PPA test.
- 2. Negative:** No colour change (medium remains straw/yellow color; no PPA to react with ferric chloride).

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
*Pseudomonas* (+) od *Enterobacteriaceae* (-)


The oxidase test is designed for specifically detecting the presence of the terminal enzyme system in aerobic respiration called cytochrome C oxidase or cytochrome a<sub>3</sub>. Cytochrome C oxidase is the terminal or last H<sub>2</sub> electron acceptor in aerobic respiratory mechanism which is composed of a number of enzymes which alternatively oxidize and reduce each other by donating or accepting electrons derived from H<sub>2</sub>.

## Koagulasa

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*S. aureus* (+) od *S. epidermidis* (-)

Coagulase test is used to differentiate *Staphylococcus aureus* (positive) which produce the enzyme coagulase, from *S. epidermidis* and *S. saprophyticus* (negative) which do not produce coagulase.

i.e Coagulase Negative *Staphylococcus* (CONS).

Principle of Coagulase Test

Coagulase is an enzyme-like protein and causes plasma to clot by converting fibrinogen to fibrin. *Staphylococcus aureus* produces two forms of coagulase: bound and free.

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
The enzyme catalase mediates the breakdown of hydrogen peroxide into oxygen and water. The presence of the enzyme in a bacterial isolate is evident when a small inoculum is introduced into hydrogen peroxide, and the rapid elaboration of oxygen bubbles occurs. The lack of catalase is evident by a lack of or weak bubble production. The culture should not be more than 24 hours old.

Bacteria thereby protect themselves from the lethal effect of Hydrogen peroxide which is accumulated as an end product of aerobic carbohydrate metabolism.

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


FIG. 1. Bacteria tagged with FLVPs. (A) Marine heterotrophic bacterial strain PWH3a tagged with phage isolate PWH3a-P1, which was fluorescently labeled with POPO-1 (blue fluorescence). (B) A natural marine bacterial community which includes cyanobacteria of the genus *Synechococcus* (red autofluorescence) as well as marine *Synechococcus* sp. isolate BBC1. *Synechococcus* isolate BBC1 has been tagged with the POPO-1-labeled cyanophage S-BIS1 (blue fluorescent halo). (C) A natural marine bacterial community to which *Synechococcus* sp. isolate BBC1 and a heterotrophic bacterium (isolate PWH3a) have been added. BBC1 has been tagged with the cyanophage S-BIS1 (blue fluorescent halo), and PWH3a has been tagged with phage isolate PWH3a-P1, which was fluorescently labeled with YOYO-1 (green fluorescence). The samples were viewed by epifluorescence microscopy with a blue (A) or violet (B and C) filter set, as described in the text. Pictures were taken at  $\times 1,000$  magnification (10- to 30-s exposures; Kodachrome 400 film). Scale bars, 10  $\mu\text{m}$ .

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

**Fig. 2.** T7-bio phage bound to streptavidin-functionalized QDs. TEM images of phage or phage-QD targeted bacteria are shown. The arrowheads point to QDs conjugated to the phage head. (*Inset*) Control T7-myc phage that are not biotinylated and therefore have no conjugated QDs. (Scale bars: 50 nm.)

**Fig. 4.** Fluorescence microscope images of phage-QD complexes bound to cells. (*a* and *b*) A typical picture of *E. coli* cells exposed at low multiplicity to QD-tagged biotinylated phage. The field was simultaneously illuminated with a low-intensity white light source and a fluorescence excitation ( $447 \pm 15$  nm). The images are of two different quantized blinking states of a single QD: off (*a*) and on (*b*). (*c*) Fluorescence micrograph of cells with 100-fold excess of biotinylated phage. (*d*) Bright-field transmission micrograph of the same sample area obtained immediately after capturing the image in *c*. Note that some cells are immobilized on a substrate, but some are mobile in solution, resulting in out-of-focus fluorescence images when the focus is maintained on the cells on a substrate surface. [Scale bars:  $1 \mu\text{m}$  (*a* and *b*) and  $2 \mu\text{m}$  (*c* and *d*).]

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TABLE I. Properties of the nucleic acid stains.<sup>32,33</sup>

Stain	$\lambda_{Ex}$	$\lambda_{Em}$	Enhancement factor <sup>a</sup>	Quantum yield	Binding properties <sup>b</sup>
DAPI	358	461	20	—	Semi-permeant; AT selective; binds to dsDNA
SYBR gold	300, 495	537	>1000	0.6–0.7	Permeant; binds to RNA, ssDNA, and dsDNA
Ethidium bromide	518	605	20–30	0.15	Impermeant; binds to RNA, ssRNA, sdDNA, trDNA
YOYO-1	491	509	100–1000	0.52	Impermeant; binds to ssDNA and dsDNA

<sup>a</sup> Fluorescence enhancement of stain after binding to the nucleic acid.

<sup>b</sup> ss = single stranded, ds = double stranded, ts = triple stranded; AT = adenine-thymine.


FIG. 2. SYBR gold image of *Salmonella typhimurium* LT2 (rod shapes) and phage P22 (bright spots). Magnification = 1000 $\times$ .

**Table 1** Phage-based assays for the detection of foodborne pathogens

Targeted pathogen	Phage (host range)	Sample matrix	Detection scheme	Detection limit	Response time
<i>Escherichia coli</i>	$\lambda$ (broad)	keberg lettuce	Reporter gene ( <i>lacI</i> and <i>lacR</i> )	1 CFU mL <sup>-1</sup>	10.3 h
	$\lambda$ Charon 30 (broad)	Sterile milk	Reporter gene ( <i>lacAB</i> )	10 cells mL <sup>-1</sup>	1.5 h
<i>E. coli</i> O157:H7	CSLO157 (narrow)	Ground beef	Phage-mediated cell lysis	<10 cells 25 g <sup>-1</sup>	10 h
	LG1 (broad)	Ground beef	Labeling of phage DNA (YOYO-1 fluorescent dye)	2 CFU g <sup>-1</sup>	6 h
	LG1 (broad)	Raw milk	Labeling of phage DNA (YOYO-1 fluorescent dye)	10 CFU mL <sup>-1</sup>	10 h
	LG1 (broad)	Ground beef	Phage amplification	2 CFU 25 g <sup>-1</sup>	23 h
	PP01 (narrow)	Apple juice	Reporter gene ( <i>lacI</i> and <i>lacR</i> )	1 CFU mL <sup>-1</sup>	22 h
	PP01 (narrow)	Spinach rinsate	Reporter gene ( <i>lacI</i> and <i>lacR</i> )	1 CFU mL <sup>-1</sup>	6 h
	T4 (broad)	Beef sirloin tip	Reporter gene ( <i>lacZ</i> , colorimetric substrate)	10 <sup>1</sup> CFU 100 cm <sup>-2</sup>	12 h
	T4 (broad)	Beef sirloin tip	Reporter gene ( <i>lacZ</i> , luminescent substrate)	100 CFU 100 cm <sup>-2</sup>	10 h
Enteric indicator bacteria	Cocktail of 3 phages, including $\lambda$ Charon 30 (broad)	Swine carcasses and slaughterhouse surfaces	Reporter gene ( <i>lacAB</i> )	10 cells cm <sup>-2</sup> or g <sup>-1</sup>	5 h
<i>Listeria monocytogenes</i>	AS11 (broad)	Shrimp, milk, cottage cheese, cabbage, lettuce	Reporter gene ( <i>lacAB</i> )	1 CFU g <sup>-1</sup>	20 h
	AS11 (broad)	Liverwurst, soft cheese	Reporter gene ( <i>lacAB</i> )	10 CFU g <sup>-1</sup>	20 h
	AS11 (broad)	Hard cheese, ground meat	Reporter gene ( <i>lacAB</i> )	1-10 CFU g <sup>-1</sup>	44 h
	AS11 (broad)	Chocolate pudding, Ricotta cheese	Reporter gene ( <i>lacAB</i> )	0.1 CFU g <sup>-1</sup>	20 h
	CBD-118 (narrow)	Turkey breast and ground meat	Phage components (cell wall binding domains)	1 CFU g <sup>-1</sup>	6 h
	CBD-118 (narrow)	Salmon, cheese, iceberg lettuce, milk	Phage components (cell wall binding domains)	10 CFU g <sup>-1</sup>	6 h
	CBD-118 (narrow)	Raw milk	Phage components (cell wall binding domains)	10 <sup>2</sup> CFU mL <sup>-1</sup>	50 h
	CBD-500 (narrow)	Soft cheese	Phage components (cell wall binding domains)	1 CFU g <sup>-1</sup>	6 h
	CBD-500 (narrow)	keberg lettuce, cheese, salmon, milk	Phage components (cell wall binding domains)	10 CFU g <sup>-1</sup>	6 h
	CBD-500 (narrow)	Turkey breast, ground meat	Phage components (cell wall binding domains)	100 CFU g <sup>-1</sup>	6 h
<i>Salmonella</i>	Felix-01	Milk	Phage amplification	<5 cells mL <sup>-1</sup>	24 h
	P22 (narrow)	Chicken breast	Phage amplification	10 <sup>4</sup> CFU g <sup>-1</sup>	3-4 h
	P22 (narrow)	Milk and eggs	Reporter gene ( <i>inuH</i> )	10 cells mL <sup>-1</sup>	2 h
	SJ2 (narrow)	Milk powder, chicken rinses, ground beef	Phage amplification	3 CFU 25 g <sup>-1</sup> or 25 mL <sup>-1</sup>	20 h

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**Table 2** Phage-based assays for the detection of pathogens in processed waters and wastewaters

Targeted pathogen	Phage (host range)	Sample matrix	Detection scheme	Detection limit	Response time	References
<i>E. coli</i>	IP008e-/2xGFP (narrow)	Sewage	Reporter gene (GFP)	NR	6 h	[17]
	IP052e-/2xGFP (narrow)	Sewage	Reporter gene (GFP)	NR	6 h	[17]
	T4 (broad)	Canal water	Labeling of phage DNA (DAPI fluorescent dye)	NR	30 min	[8, 113]
	T4 (broad)	Mineral water	Biosensor (electrochemical impedance sensing)	$10^4$ CFU mL <sup>-1</sup>	1 h	[81, 113]
	T4e- (narrow)	Sewage	Reporter gene (GFP)	NR	1 h	[16, 79]
<i>E. coli</i> O157:H7	PhiX174 (broad)	River water	Phage-mediated cell lysis	$6 \times 10^5$ CFU mL <sup>-1</sup>	2 h	[80, 114]
	PP01 (narrow)	Tap water	Reporter gene ( <i>luxI</i> and <i>luxR</i> )	1 CFU mL <sup>-1</sup>	12.5 h	[63]
	<i>Micrococcus phosphovorius</i>	ΦMP1 (narrow)	Activated sludge	Labeling of phage DNA (SYBR Green fluorescent dye)	$\sim 10^2$ cells mL <sup>-1</sup>	25 min
ΦMP2 (narrow)		Activated sludge	Labeling of phage DNA (SYBR Green fluorescent dye)	$\sim 10^2$ cells mL <sup>-1</sup>	25 min	[82]
<i>Bacillus anthracis</i>	JRB7 (narrow)	Distilled water	Phage-display biosensor (magnetoelastic)	$10^3$ CFU mL <sup>-1</sup>	20 min	[83]
	JRB7 (narrow)	Distilled water	Phage immobilization	$10^4$ spores mL <sup>-1</sup>	1 h	[85]
<i>B. anthracis</i> and <i>Salmonella</i> Typhimurium (simultaneous detection)	JRB7 (for <i>B. anthracis</i> ) and E2 (for <i>Salmonella</i> Typhimurium) (narrow)	Distilled water	Phage-display biosensor (magnetoelastic)	$5 \times 10^3$ CFU mL <sup>-1</sup>	$\sim 1$ min	[84]

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**Table 3** Phage-based assays for epidemiological and clinical diagnostics

Targeted pathogen	Phage (host range)	Sample matrix	Method	Detection limit	Response time
<i>Mycobacterium tuberculosis</i>	phAE142	Sputum	Reporter gene ( <i>luc</i> )	$10^3$ CFU mL <sup>-1</sup>	1-2 weeks
	D29	Sputum	Phage amplification	100-300 cells mL <sup>-1</sup>	48 h
	(FASTPlaque™) (narrow) D29 (FASTPlaque-Response™) (narrow)	Sputum	Phage amplification	Antibiotic resistance profiling	48 h
<i>Yersinia pestis</i>	Proprietary (KeyPath™)	Blood	Phage amplification	Antibiotic resistance profiling	5.5 h
	ΦA1122	Blood	Phage amplification	$10^6$ CFU mL <sup>-1</sup>	5 h
<i>Staphylococcus aureus</i>	ΦA1122	Serum	Reporter gene ( <i>luxAB</i> )	$10^6$ CFU mL <sup>-1</sup>	<30 min
	Proprietary	Blood	Phage amplification	10 CFU mL <sup>-1</sup>	5.5 h
	Phage-displayed peptide	Blood	Phage display and quantum dots	$10^3$ CFU mL <sup>-1</sup>	NR
<i>E. coli</i>	λ Charon 30	Urine	Reporter gene ( <i>lux</i> )	10-100 cells mL <sup>-1</sup>	1 h