Surface Analysis - XPS

ANTON MANAKHOV PHD IN PHYSICS




Central European Institute of Technology BRNO | CZECH REPUBLIC



MASARYK UNIVERSITY Czech Republic

Why we need to know surface chemistry?

- Surface properties affects interaction of material (metals, plastic, powder, etc.) with the environment (water, polymers, in vivo...) and specific chemical composition and morphology is required for particular application.
- Surface Chemistry influences cell adhesion, fouling of bioorganisms, adhesion between different phases of composite materials.
- Often the application of material with noble bulk properties (hardness, flexibility, chemical stability) is not possible due to low surface free energy, or too reactive surface instead.
- Deposition of thin film can tune surface properties without degradation of the bulk of the material
 Roughness
 Electronic properties



Surface analysis methods :

| Technique | Probe In/Out | Depth resolu- tion, nm | Lateral resolu- tion, nm | Informa- tion | Advantage (+) and drawback (-) |
|---|----------------------------|---------------------------------|--------------------------------|---|---|
| Contact angle | Liquid water droplet | 0.1 | 1000 | Surface energy | + Fast acquisition - No molecular information |
| X-Ray photoelectron spectroscopy XPS | X-Ray/ electrons | 5 | 3000 | Elemental composition (except H), binding state | + Quantitative + Information about the neighbours of the atoms - Limited Molecular information - sensitivity >0.1 at% |
| Scanning electron microscopy , SEM | Electrons/ electrons | 3 | 2 | Surface image | - No direct information regarding the topography. |

| Technique | Probe In/Out | Depth resolu- tion, nm | Lateral resolu- tion, nm | Informa- tion | Advantage (+) and drawback (-) |
|---|---------------------|---------------------------------|--------------------------------|--|---|
| Energy Dispersive analysis | Electrons/ X-Ray | 2000 | 100 | Elemental information | - Low precision of quantification |
| IR attenuated total reflection ATR-FTIR | IR/IR | ~2000 | 2000 | Surface composition binding state | + Fast acquisition +Information regarding the molecular functions - Overlap of some chemical domains -Quantification is not possible |
| Time of Flight Secondary Ion Mass Spectrometry (ToF-SIMS) | Ions/Ions | 1 | 100 | Surface composition | + Highly Sensitive (>1 ppm) +Molecular information is possible + 3D reconstruction of the layer - non quantitative |


Methods based on the radiation of material by photons : XPS

- XRD
- UPS
- UV-Vis spectrophotomentry
 - IR spectroscopy

Interaction of photons with matter X-ray X-ray The interaction ionization X-rays ionization of radiation with Compton matter. Scattering Longer wavelength Photoionization X-ray Ionization Ultraviolet Electron energy •^^^^ level changes. Large number of available energy Visible states, strongly absorbed. Infrared Molecular vibration Molecular Small number of Microwaves rotation available states. nd torsion almost transparent.


X-ray Photoelectron Spectroscopy XPS

X-ray Photoelectron Spectroscopy (XPS), also known as Electron Spectroscopy for Chemical Analysis (ESCA) is a widely used technique to investigate the chemical composition of surfaces. X-ray Photoelectron spectroscopy, based on the photoelectric effect, was developed in the mid-1960's by Kai Siegbahn and his research group at the University of Uppsala






The Photoelectric Process




XPS involves irradiating a sample with X-rays of a characteristic energy and measuring the flux of electrons leaving the surface. The energy spectrum for the ejected electrons is a combination of an overall trend due to transmission characteristics of the spectrometer, energy loss processes within the sample and resonance structures that derive form electronic states of the material under analysis. The instrumental contribution is an unwelcome fact of the measurement process, but the background and resonance peaks offer information about the sample surface

Auger Relation of Core Hole




Emitted Auger Electron


- L electron falls to fill core level vacancy (step 1).
- KLL Auger electron emitted to conserve energy released in step 1.
- The kinetic energy of the emitted Auger electron is: KE=E(K)-E(L2)-E(L3).



 $h \cdot v = E_b + E_{kin} + \Phi$




Because the Fermi levels of the sample and spectrometer are aligned, we only need to know the spectrometer work function, Φ_{spec} , to calculate BE(1s).




Chemical Shifts

Chemical Shifts- Electronegativity Effects

| Functional Group | | Binding Energy (eV) |
|---------------------|------------------------------|------------------------|
| hydrocarbon | <u>С</u> -Н, <u>С</u> -С | 285.0 |
| amine | <u>C</u> -N | 286.0 |
| alcohol, ether | <u>С</u> -О-Н, <u>С</u> -О-С | 286.5 |
| CI bound to C | <u>C</u> -Cl | 286.5 |
| F bound to C | <u>C</u> -F | 287.8 |
| carbonyl | <u>C</u> =O | 288.0 |





Heavy elements

 $L_{1} + L_{2} = L$ $\vec{2} + \vec{2} = 0, 1, 2, 3, 4$ $\vec{S}_{1} + \vec{S}_{2} = \vec{S}$ $\vec{1/2} + \vec{1/2} = 0, 1$ $\vec{J} = \vec{L} + \vec{S} = 0, 1, 2, 3, 4$


Final Effects – shake-up

Final State Effects - Shake-up/ Shake-off

Results from energy made available in the relaxation of the final state configuration (due to a loss of the screening effect of the core level electron which underwent photoemission).

Shake-up: Relaxation energy used to excite electrons in valence levels to bound states (monopole excitation).


Shake-off: Relaxation energy used to excite electrons in valence levels to unbound states (monopole ionization).



Electron Scattering Effects Energy Loss Peaks


$$e_{ph} + e_{solid} \rightarrow e_{ph}^{*} + e_{solid}^{**}$$

Electron Scattering Effects Plasmon Loss Peak





In metals the free electrons are constrained to move within energy bands that are characteristic of the material and these material properties influence the shape of the energy loss distribution, namely, scattering of the photoelectric electrons by free electrons with discrete energy bands produces energy loss distributions with relatively narrow structures.


Depth of analysis and Mean Free Path



from www.lasurface.com










Quantification


$$I_{\rm A}^i \sim Q \mathcal{A} \ c_{\rm A} \sigma_{\rm A}^i \lambda^i T^i L_{\rm A}^i f(\phi, \theta), \tag{7.25}$$


where Q is the photon flux $[cm^{-1}s^{-1}]$, \mathcal{A} effective are of the sample, c_A concentration of A, σ_A^i is the partial ionization cross section, λ^i mean free path, T^i the transmission, L_A^i angular asymmetry coefficient and $f(\phi, \theta)$ function, depending on the geometry of the experiment.

















© 2010 WILEY-VCH Verlag GmbH & Co. KGaA, Weinheim



AllylAmine Plasma Polymers






High concentration of amines in AA pp can be obtained only at low plasma power but it leads to huge loss of material.

Effects of amide and amine plasma-treated ePTFE vascular grafts on endothelial cell lining in an artificial circulatory system

David Y. Tseng,¹ Elazer R. Edelman^{1,2}

¹Harvard–MIT Division of Health Sciences and Technology, Massachusetts Institute of Technology, Cambridge, Massachusetts 02139 ²Department of Medicine, Cardiovascular Division, Brigham and Women's Hopsital, Harvard Medical School, Boston, Massachusetts 02115






Figure 6. Bar graph showing the percentages of endothelial cell adhesion on plasma coated (open bar) and uncoated (solid bar) PTFE disks, respectively. Statistical analysis showed a significant difference of endothelial cell adhesion between uncoated and coated PTFE disks (p < 0.005, n = 8). The data are presented as mean \pm SE.

