

No.	Name	Type	Ct
1	10x UBC10	Standard	14.88
2	10x UBC10	Standard	14.82
3	100x UBC10	Standard	17.54
4	100x UBC10	Standard	17.6
5	1 000x UBC10	Standard	22.32
6	1 000x UBC10	Standard	22.37
7	10 000x UBC10	Standard	25.41
8	10 000x UBC10	Standard	26.23
9	100 000x UBC10	Standard	29.04
10	100 000x UBC10	Standard	28.93
11	10x LBD29	Standard	18.95
12	10x LBD29	Standard	19.71
13	100x LBD29	Standard	22.87
14	100x LBD29	Standard	22.54
15	1000x LBD29	Standard	27.23
16	1000x LBD29	Standard	26.57
17	10 000x LBD29	Standard	29.35
18	10 000x LBD29	Standard	
19	100 000x LBD29	Standard	
20	100 000x LBD29	Standard	36.58
21	DMSO 6H UBC10	Unknown	14.88
22	DMSO 6H UBC10	Unknown	13.96
23	DMSO 6H UBC10	Unknown	15.06
24	NAA 6H UBC10	Unknown	14.78
25	NAA 6H UBC10	Unknown	14.9
26	NAA 6H UBC10	Unknown	14.89
27	DMSO 12H UBC10	Unknown	15.41
28	DMSO 12H UBC10	Unknown	14.89
29	DMSO 12H UBC10	Unknown	13.91
30	NAA 12H UBC10	Unknown	14.86
31	NAA 12H UBC10	Unknown	14.82
32	NAA 12H UBC10	Unknown	15.93
33	DMSO 6H LBD29	Unknown	23.81
34	DMSO 6H LBD29	Unknown	23.67
35	DMSO 6H LBD29	Unknown	23.87
36	NAA 6H LBD29	Unknown	18.73
37	NAA 6H LBD29	Unknown	18.69
38	NAA 6H LBD29	Unknown	18.8
39	DMSO 12H LBD29	Unknown	24.49
40	DMSO 12H LBD29	Unknown	22.97
41	DMSO 12H LBD29	Unknown	23.87
42	NAA 12H LBD29	Unknown	19.72
43	NAA 12H LBD29	Unknown	19.77
44	NAA 12H LBD29	Unknown	19.8
45	DMSO 24H UBC10	Unknown	16.66
46	DMSO 24H UBC10	Unknown	16.62
47	DMSO 24H UBC10	Unknown	16.59
48	NAA 24H UBC10	Unknown	15.7
49	NAA 24H UBC10	Unknown	15.78

\*Ct na úrovni negatívnej kontroly (preto chýb)

\* Ct na úrovni negatívnej kontroly (preto chýb)

50	NAA 24H UBC10	Unknown	15.79
51	DMSO 48H UBC10	Unknown	15.8
52	DMSO 48H UBC10	Unknown	16.71
53	DMSO 48H UBC10	Unknown	16.6
54	NAA 48H UBC10	Unknown	17.57
55	NAA 48H UBC10	Unknown	17.55
56	NAA 48H UBC10	Unknown	17.52
57	DMSO 24H LBD29	Unknown	25.48
58	DMSO 24H LBD29	Unknown	24.52
59	DMSO 24H LBD29	Unknown	25.51
60	NAA 24H LBD29	Unknown	18.8
61	NAA 24H LBD29	Unknown	18.65
62	NAA 24H LBD29	Unknown	18.87
63	DMSO 48H LBD29	Unknown	24.67
64	DMSO 48H LBD29	Unknown	24.88
65	DMSO 48H LBD29	Unknown	23.93
66	NAA 48H LBD29	Unknown	20.83
67	NAA 48H LBD29	Unknown	19.85
68	NAA 48H LBD29	Unknown	19.96
69	NEG UBC10	NTC	39.16
70	NEG LBD29	NTC	
71	NEG UBC10	NTC	
72	NEG LBD29	NTC	39.17

a hodnota) - špatne napipetovaný templát? - pri riedení 10000x použite teda iba hodnotu zo vzorky č 17  
a hodnota) - dáta z riedenia 100 000x do konštrukcie kalibračnej krivky LBD29 vynechajte



